Estimado cliente

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Dear Customer

Thank you for choosing this high-quality product from Dentaureum.

In order to achieve the best results with this product, it is important to carefully study and follow these instructions for use.

The written instructions cannot cover all possible situations that can occur during use. For this reason, our hotline is available to answer any other questions or concerns that may arise.

Due to constant developments, we recommend, even when using the same products regularly, that you study the enclosed current instructions for use or refer to our website at www.dentaureum.de.

1. Manufacturer

Dentaureum GmbH & Co. KG
Turnstr. 31 | 75228 Ispringen | Germany

2. General product description

The hyrax® click screw (hygienic rapid expansion screw) is a special expansion screw with a patented locking device (Fig. 1) that prevents unwanted turning back during quick mid-palatal suture separation. The retention legs of the palatal expansion screw are welded or soldered onto preformed metal bands (Fig. 2) or cap splints (Fig. 3). Alternatively, they can be polymerized into acrylic splints (Fig. 4). The bands or splints are cemented in the usual way.

3. Area of application

Used for rapid palatal expansion with a fixed appliance.

4. Fabricating a hyrax® click appliance using the example of preformed bands

Preformed bands are fitted on the maxillary left and right 1st bicuspsids and left and right 1st molars. These bands are transferred to a dental stone model using an impression or fitted to a dental stone model. The legs of the hyrax® click screw are then bent to fit the bands so that they can be soldered or welded into place (Fig. 2 + 5).

The hyrax® click screw should "float" freely at least 1-2 mm above the palate roof/gingiva (Fig. 6). Its dorsal side should form a straight line with the mesio palatal cusps of the 1st molars. To avoid any leverage effect, be sure that the occlusal level is parallel. The hyrax® click screw is positioned so that it may be turned in a dorsal direction (Fig. 2 – 4) in order to make the task of inserting the safety key and activating the screw easier.

5. Information on use of hyrax® click screws

To increase strength, the bands can be joined on the buccal and/or lingual side with a piece of 1.0 mm remalay® wire (REF 528-100-00). Laser weld or solder these wires above the legs. The palatal strengthening wire should project slightly beyond the bands and be supported by the adjacent teeth to allow better distribution of expansion forces.

In order to bend the retention arms of the hyrax® click screw gently and easily, we recommend using the universal bending device (REF 611-122-00, Fig. 7 + 9) without the aid of pliers (lever principle). If this is not available, the retention arms must be stabilised by holding at the expansion screw body with a pair of sturdy pliers e.g. Wire bending pliers Maxi (REF 013-522-00), in order to prevent the retention arms from breaking (Fig. 8).

6. Composition

The composition is included in the materials list, refer to www.dentaureum.de.

7. Safety information

Always use the blue safety key (Fig. 10) supplied to prevent swallowing!

8. Information on single-use products

The product is intended for single use only. Reconditioning of a product that has already been used (recycling) and its reuse on a patient is not permitted.

9. Availability

Description	straight legs	bent legs	Quantity
hyrax® click Mini - 7 (Expansion 7 mm)	602-830-10 602-830-30	602-831-10 602-831-30	1 piece 10 pieces
hyrax® click Medium - 10 (Expansion 10 mm)	602-832-10 602-832-30	602-833-10 602-833-30	1 piece 10 pieces
hyrax® click Maxi - 12 (Expansion 12 mm)	602-834-10 602-834-30	602-835-10 602-835-30	1 piece 10 pieces

with safety key

Safety key	611-120-01	10 pieces
Universal bending tool	611-122-00	1 piece

10. With reference to our quality standards

Dentaureum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

11. Explanation of symbols used on the label

⚠ Please refer to the label. Additional information can be found at www.dentaureum.de (Explanation of the symbols used on the label REF 989-313-00).

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